


# GeneProof®

## GPack

### Covid-19 • Flu • RSV

Differential diagnostics of many pathogens from one sample in a single run



**Simplify your laboratory workflow!**

One sample  
One extraction  
One PCR run



many  
results

**NEW  
GPack**

#### **Covid-19 • Flu • RSV GPack contains:**

- GeneProof SARS-CoV-2 Advanced PCR Kit, 100 rxn
- GeneProof Flu Multiplex PCR Kit, 100 rxn

#### **1 SAMPLE - 5 EXAMINED PATHOGENS**

- Diagnostics of SARS-CoV-2, *Influenza A*, *Influenza B*, RSV A/B
- The same symptoms - different diagnosis

#### **EXCELLENT QUALITY**

- Quality proven by EQA panel results
- Unique specificity for SARS-CoV-2

#### **UNIVERSAL PCR PROFILE**

- Enables the combinability of PCR kits from different diagnostic groups in one run
- Simplifies laboratory workflow and saves time

#### **TECHNOLOGICALLY ADVANCED PRODUCTS**

##### **Easy-to-Use concept**

- Single tube Ready-to-Use Master Mix
- No additional PCR reagents pipetting necessary

##### **Contamination prevention**

- dUTP-UNG contamination prevention technology

##### **False negative results control**

- Unique internal control for the control of the whole diagnostic process including RNA extraction, reverse transcription and PCR amplification
- Endogenous internal control for the control of the proper sampling



COMPATIBLE WITH A WIDE  
RANGE OF REAL-TIME PCR  
DEVICES



CERTIFIED  
DIAGNOSTIC TESTS